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REV GA-03-SIBD-XX/2021-05

Besmed Health Business Corp.  
No. 5, Lane 116, Wu-Kang 2nd Road,  
Wu-Ku District  
24888 New Taipei City  
Taiwan (R.O.C.)

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<b>1</b>	General .....	Page	<b>11</b>
<b>2</b>	Safety and Responsibilities .....	Page	<b>11</b>
<b>3</b>	Modification of the product .....	Page	<b>12</b>
<b>4</b>	Before using the product .....	Page	<b>12</b>
<b>5</b>	Using the product .....	Page	<b>13</b>
<b>6</b>	Product Maintenance / Reprocessing .....	Page	<b>13</b>
<b>6.1</b>	Warnings .....	Page	<b>13</b>
<b>6.2</b>	Restrictions due to reprocessing .....	Page	<b>13</b>
<b>6.3</b>	Preparation at the place of use .....	Page	<b>14</b>
<b>6.4</b>	Preparation for cleaning .....	Page	<b>14</b>
<b>6.5</b>	Manual cleaning and disinfection .....	Page	<b>14</b>
<b>6.6</b>	Automatic cleaning and disinfection .....	Page	<b>14</b>
<b>6.7</b>	Drying .....	Page	<b>15</b>
<b>6.8</b>	Maintenance, control and inspection .....	Page	<b>15</b>
<b>6.9</b>	Packaging .....	Page	<b>15</b>
<b>6.10</b>	Sterilisation .....	Page	<b>15</b>
<b>6.11</b>	Additional Information .....	Page	<b>16</b>
<b>7</b>	Accessories .....	Page	<b>16</b>
<b>8</b>	Replacement parts and repairs .....	Page	<b>16</b>
<b>9</b>	Disposal .....	Page	<b>16</b>
<b>10</b>	Explanation of symbols .....	Page	<b>17</b>

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## 1 General

- This product may only be used by trained professionals under medical instruction.
- The product is intended for manual monitoring and ventilation during artificial ventilation.
- It is mainly used in hospitals in the operating room for anaesthesia and machine ventilation.
- These instructions reflect the state of the art at the time of publication.
- This product complies with the provisions of the EU Directive 93/42/EEC (Medical Device Directive) and Medical Device Regulation (EU) 2017/745.
- The instructions for use comply with DIN EN 82079-1.

## 2 Safety and Responsibilities

**IMPORTANT! READ CAREFULLY BEFORE USE**  
**KEEP FOR FUTURE REFERENCE**

- Please observe the following instructions and care instructions before use.
- If in doubt about the use and connection of the product please contact customer service.
- The user should have sufficient experience and knowledge in the assembly, disassembly, cleaning and the use of the product.
- If additional oxygen is administered care should be taken to ensure that smoking and other naked flames are forbidden in the vicinity of oxygen-carrying equipment and fittings.
- The product must not come into contact with a pointed or a sharp object which could damage it. This could lead to breakages or similar damage and the product would start to leak.
- Keep the produce free of oil and grease, otherwise there is a risk of explosion. This includes the exposure to cleaning alcohols, hand cream and adhesives – in contact with highly compressed oxygen explosive reactions may occur.
- The product may not be used for ventilating babies.
- Make sure that the correct size is used for each individual patient.
- Reusable products are delivered non-sterile and must therefore be unpacked, disinfected and sterilised (according to chapter 6) before use.
- The processing instructions under chapter 6 must be followed carefully to avoid any exposure to danger! In case of questions please contact customer service.
- Repeated use can lead to hardening and cracking. In this case the product must be disposed of according to chapter 9.

- When used in patients with a suspected prion disease, there may be a high risk of infection.

It lies within the discretion of the physician to properly dispose of the product according to chapter 9 or to reprocess it according to national regulations (see chapter 6).

### 3 Modification of the product

- Apart from the described maintenance instructions, the product may not be modified or repaired by the operator!
- The manufacturer's responsibility is limited to the dimensional and physiological suitability of the product.
- Improper use, modifications to the product or negligence will void any warranty. This also applies to unauthorised repairs. Normal signs of wear and tear on the device are also excluded.

### 4 Before using the product

- Do not expose the product to direct sunlight or other light sources. They must be stored and transported in a dry place and in their original packaging.
- Remove all delivered parts from the packaging.  
If necessary, remove the respective plastic covers and dispose of the packaging according to the current national regulations.
- Before each use verify that parts or components do not show any signs of hardening, damage, breakage or cracks. The product must also not be contaminated or soiled. If any of the conditions described above are present, the device must not be put into operation and must be disposed of in accordance with chapter 9.
- In case of doubt perform a comparison check with a brand-new product.  
**→ Damaged products must not be used!**  
They must be disposed of in accordance with the instructions in chapter 9.
- The product contains silicone.
- The product is delivered in a non-sterile condition. Cleaning and disinfection must be performed before each use according to instructions in chapter 6.
- To attach and detach to/from the system, the product must be held at the connector to prevent cracks from forming between the connector and the main tubing.
- Attach the product to a suitable connection tube.
- After assembly, the product must be checked for sufficient pull-off force and for kinks as well as any cracks or hardening.

## 5 Using the product

- After connection to a respiratory circuit, the product is ready for use.
- Choose the correct type and size of mask according to the anatomical features of the patient.  
If applicable observe national legislation.
- Fill the breathing bag with oxygen.
- Compress the bag according to medical instruction.

## 6 Product maintenance/Reprocessing

- The product is reusable subject to risk assessment by the qualified user.
- However, the service life is limited to a maximum of **two years or 100 reprocessing cycles** from the first usage onwards.  
The actual wear and tear and the associated number of reprocessing cycles will vary. For this reason, a functional test must be carried out before each use to determine possible material fatigue or impaired functionality.
- Observe the following procedures for reprocessing (cleaning and disinfecting). These are validated by the manufacturer.
- A complete reprocessing cycle always consists of cleaning and disinfection (manual or mechanical) and sterilization.

### 6.1 Warnings

- The product is delivered non-sterile.
- The product must be processed before initial use and before each repeat use.  
→ **The delivered packaging cannot be sterilized!**
- The products may only be processed by specifically trained personnel.
- Do not use any pointed or sharp objects as well as hard brushes, scouring agents or other materials that could damage the surface materials!
- The manufacturer accepts no responsibility for damage due to incorrect processing.
- The processing of the product lies within the clinical responsibility of the user.
- No residues may remain in or on the product.
- Do not use damaged products.  
If necessary, contact customer service.

### 6.2 Restrictions due to reprocessing

- The use of oils and greases is not recommended for silicone products.
- Do not use disinfectants based on formaldehyde, glutaraldehyde, glutoprotamine, coco dipropylene triamine, quaternary ammonium compounds, phenol derivatives or iodine-cleaving compounds.
- In addition, it must be ensured that no hard or sharp-edged objects come into contact with the product during disinfection and that the product is rinsed thoroughly afterwards.

### 6.3 Preparation at the place of use

- Immediately after use, manual disinfection should be performed to reduce the risk of infection for users.
- It is recommended to rinse the product thoroughly with warm water and a diluted sodium bicarbonate solution (8-10% by volume) immediately after use until all visible foreign matter has been removed.
- Do not let any residues become dried-on.

### 6.4 Preparation for cleaning

- Remove the product from the connecting device.

### 6.5 Manual cleaning and disinfection

- Solely manual reprocessing is not recommended due to the lower effectiveness and reproducibility.

It should only be used if no automated procedure is available.

- Use fully demineralised water.
- Cleaning and disinfecting solution must be prepared before each manual cleaning cycle. Please follow the manufacturer's instructions for the solution.
- A solution of Sekusept Aktiv is recommended.
- Place the product in the cleaning and disinfecting solution. Pay attention to the manufacturer's recommended concentration and soaking time.
- Even when cleaning manually, make sure that all cavities are properly rinsed through.
- A small brush with soft bristles or a soft sponge can be used to help clean the product.
- After cleaning, rinse sufficiently with fully demineralised water to remove any residues.
- Avoid any water accumulation.
- Dry the product completely.
- Check visually for residues.  
Repeat the rinsing process if necessary.

### 6.6 Machine cleaning and disinfection

- Use fully demineralised water.
- If using alkaline cleaning agents, neutralisation must be carried out first.
- Do not use any desiccants.
- It is recommended to use an automated/mechanical cleaning and disinfection procedure.
- The product is suitable for machine processing/thermal disinfection ( $\leq 93^{\circ}\text{C}$ ).
- It must be ensured that the selected washer-disinfector complies with ISO 15883.
- Follow the operating instructions of the respective device model.

- Use cleaning agents that are suitable for silicone and thermostable plastics.
- Recommended cleaning agents:  
Cleaner: Sekumatic ProClean // Neutraliser: ECOLAB Sekumatic FNZ
- The operating cycle (cleaning, disinfection, rinsing, drying if necessary) must achieve the specified performance.
- Do not overload the basket and avoid shadowing.
- Load the product with openings and recesses with open side facing downwards to avoid accumulation of cleaning water.

#### **Instructions:**

1. Place the product in the washer-disinfector / arrange around the spiral rack (ensure that all cavities are effectively flushed through).  
Observe the loading pattern of your washer-disinfector if necessary.
2. Start a cycle with the following parameters: (e.g. VARIO TD)
 

– Cleaning with alkaline cleaner	max. 55°C // 5–10 min.
– Neutralisation	20°C // 2–6 min.
– Intermediate rinsing	20°C // 2–3 min.
– Disinfection	93°C // 5 min.
– Drying	max. 100°C // 30 min.
3. Inspect for any residual soil. Repeat the cycle if necessary.
4. Inspect the product and prepare for sterilization.

### **6.7 Drying**

- After cleaning, the products must be dried thoroughly (preferably supported by compressed air).
- Do not use drying agents or cloths which can leave residues on the silicone.

### **6.8 Maintenance, control and inspection**

- Inspect the products after disinfection and before each use for intactness, leakage and cleanliness.
- Check the product for cracks and material changes.
- If residual contamination is still visible, repeat the entire procedure until no more residues are visible.
- After reprocessing and before each use, a functional check must be carried out.

### **6.9 Packaging**

- After processing, package the product according to ISO 7153.
- Store the sterilised product in sterile packaging in a closed cabinet, protected from dust, light, moisture and temperature fluctuations.

### **6.10 Sterilisation**

- The product must be sterilised before each use.
- Before sterilisation, place the product in a suitable container or in suitable sterilisation wrapping.

- The correct size of container or wrapping must be chosen.
- Use a sterilisation indicator for the wrapping and note the sterilisation and expiry date on the wrapping.
- The product is only suitable for steam sterilization in a fractionated vacuum process.  
Maximum exposure time: 5 min.  
Maximum temperature: 134°C
- Follow the instructions of the autoclave manufacturer carefully.
- After steam sterilisation, the product must be allowed to cool down at room temperature.
- The manufacturer accepts no liability for risks arising from improper sterilisation (contamination).

## 6.11 Additional Information

- Frequent sterilisation can accelerate the natural ageing of the product. In addition, the service life of sterilised products depends very much on carefully performed pre-treatment (cleaning and disinfection).
- When sterilising several products in one sterilisation cycle, the maximum load of the steriliser must not be exceeded.

## 7 Accessories

- The user is responsible for the correct combination with other medical products.
- The product can be used in combination with a ventilation tube with a 22 mm ISO connection.

## 8 Replacement parts and repairs

- Repairs must only be carried out by the manufacturer or by authorised persons.
- In case of functional or material problems, contact customer service immediately.
- If the cause of the error cannot be eliminated, the product must not be put into operation. After consulting the customer service, the product must be disposed of according to chapter 9 or returned to the sales distributor.
- In the event of a return for repair or replacement, the entire reprocessing process must be carried out beforehand in order to exclude any risk to the manufacturer's employees.

We reserve the right to reject soiled or contaminated products for safety reasons.

## 9 Disposal

- If any of the above mentioned errors have occurred, the product must be disposed of.
- The product must be disposed of at the latest two years after first use or following 100 reprocessing cycles.
- Safe disposal is carried out in accordance with current national regulations.

## 10 Symbolerklärung / Explanation of symbols

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Symbol for Germany

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Symbol for Poland

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Symbol for Great Britain

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Observe the operating instructions

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Do not use any grease  
or other form of lubricant

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Medical Device

Importeur

Importer