Montage- und Gebrauchsanweisung	04	DEU
Installation and operating instructions	25	ENG
Mode d'emploi	46	FRA
Istruzioni per il montaggio e l'uso	67	ITA
Instrucciones de montaje y funcionamiento	88	ESP
Instruções de instalação e de utilização	109	POR
Montage- en gebruiksaanwijzing	130	NLD
Bruksanvisning	151	SWE
取り付けおよび使用の説明書	172	JAP

D^{med}® **Triango** 60 C
D^{med}® **Triango** 60 W
D^{med}® **Triango** 60 F

Behandlungsleuchte
Treatment light
Lumière de traitement
Luce di trattamento
Luz de tratamiento
Candeeiro de Tratamento
Behandelingslicht
Manöverlampe
処置ライト

SYMBOLE; SYMBOLS; SYMBOLES; SIMBOLI; SÍMBOLOS; SÍMBOLOS; SYMBOLEN; SYMBOLER; 記号の意味

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
IMPORTANT!
THIS OPERATING MANUAL MUST BE THOROUGHLY READ BEFORE USING THE PRODUCT!
→ STORE FOR LATER USE!

CONTENTS

1.	VARIANTS AND SCOPE OF DELIVERY	26
1.1	Triango 60 C	26
1.2	Triango 60 W	26
1.3	Triango 60 F	27
2.	SAFETY INSTRUCTIONS	28
2.1	Intended use	28
2.2	User profile	28
2.3	Safety instructions	28
2.4	Warning levels	28
3.	INSTALLATION: Triango 60 C	29
3.1	Load data.....	29
3.2	Shortening the ceiling rod	29
3.3	Ceiling bracket installation	29
3.4	Ceiling bracket installation	30
3.5	Ceiling arm installation.....	31
4.	INSTALLATION: Triango 60 W	32
4.1	Load data.....	32
4.2	Wall bracket installation	32
4.3	Wall arm installation.....	33
5.	INSTALLATION: Triango 60 F	34
6.	INSTALLATION: Triango 60 light head	36
7.	OPERATION	36
8.	CLEANING AND DISINFECTION	37
8.1	Sterilisation of the handle	38
9.	SAFETY CONTROLS	38
9.1	SETTING SPRING FORCE.....	38
10.	DEINSTALLATION	39
10.1	Disposal	39
11.	ACCESSORIES	39
12.	ADDITIONAL INFORMATION	39
13.	TROUBLESHOOTING	40
14.	TECHNICAL DATA	40
15.	NAMING THEME	42
16.	ELECTROMAGNETIC COMPATIBILITY (EMC)	43

1. VARIANTS AND SCOPE OF DELIVERY

1.1 Triango 60 C

A: Ceiling arm


B: Light head with sterilisable handle

C: Ceiling rod

D: Retaining ring

E: Ceiling plate and hood

1.2 Triango 60 W

A: Light head with sterilisable handle


B: Wall arm

C: Wall angle

D: Retaining ring

E: Wall bracket and hood

1.3 Triango 60 F

A: Light head with sterilisable handle

B: Spring arm

C: Cable holder

D: Upright rod

E: Roller stand

2. SAFETY INSTRUCTIONS

2.1 Intended use

The Triango 60 light is a small surgical light. It is a single light in the patient environment for use in treatment rooms to support diagnosis or treatment that does not pose a hazard to the patient in the event of interruption due to a light failure. It is intended for continuous operation and is not intended to be combined with other medical devices.

The essential performance of Triango 60 is the delivery of sufficient light to the operating field while limiting the delivery of radiant energy to the operating field and its observer. This is achieved by providing the capability to deliver central illuminance of at least 40 klx, while ensuring the total irradiance at the light field centre at the maximum illuminance distance does not exceed 700 W/m².

2.2 User profile

Medical specialist

All people who have completed medical training and work in their trained professional field.

Cleaning specialist

Someone who has been instructed in the national and workplace-related hygiene regulations.

Electrician

Someone who is trained in the fields of electronics and electrical engineering and knows the relevant standards and regulations.

Qualified specialist

Someone who is able to perform the installation/deinstallation due to technical training, knowledge and experience as well as knowledge of the regulations.

2.3 Safety instructions

- ▶ Operation by a medical specialist
- ▶ The manual is part of the product and must be kept and made available to all subsequent users.
- ▶ All work on the light (including repairs) may only be carried out by a qualified electrician. Installation may only be carried out by a qualified specialist.
- ▶ The light must not be changed or manipulated. Only approved original parts may be used. Use other than the intended use with the original parts can lead to different technical values and life-threatening dangers.
- ▶ Do not exceed the maximum weight, or hang from, lean on or step on it, as this may cause the device to tip over and result in serious injury.
- ▶ Operation in potentially explosive areas is prohibited. The power supply to the light represents a potential ignition source.
- ▶ The light may only be operated in dry and dust-free rooms.
- ▶ The light must not be left switched on unattended.
- ▶ Only connect the light to the mains supply with a protective conductor to avoid electric shock.

- ▶ In the case of lights of protection class I, the protective conductor must be connected to the light housing.
- ▶ A damaged light should not be used. Defective cables and a defective handle also pose a potential hazard. Do not place cables near heat sources or on sharp edges.
- ▶ Never place extra weight on the light head and the arm system.
- ▶ When in operation, the light must never be covered with a cloth or similar.
- ▶ The ventilation openings (if available) must always be free during operation!
- ▶ The light must not be operated in the vicinity of external heat sources that exceed the maximum ambient-temperature of the light.
- ▶ The light must not be used outside the intended environmental conditions.
- ▶ Do not use together with medical devices that can react sensitively to a light spectrum in the visible range (e.g. pulsating light and/or light with high illuminance)
- ▶ The light may only be used for the purpose mentioned here.
- ▶ The manufacturer cannot be held responsible for damage caused as a result of use that deviates from the intended use or the non-observance of safety instructions and warnings.
- ▶ When using several lights at the same time, the total illumination E_{Total} in the light field must not exceed 1000W/m² during operation.
- ▶ Before connecting to the mains supply, the mains data must be checked to ensure that it matches the device data.
- ▶ **Triango 60 F**
The light must be held securely during transport within the hospital.

2.4 Warning levels

DANGER

Warnings of dangers which, if the measures are ignored, can lead to **death or serious injuries**.

WARNING

Warnings of dangers which, if the measures are ignored, can lead to **injuries**.

CAUTION


Warnings of dangers which, if the measures are ignored, can lead to **damage**.

3. INSTALLATION: Triango 60 C


3.1 Load data

Bending moment M_B	135 Nm
Gravitational force F_G	140 N


3.2 Shortening the ceiling rod

- ▶ Loosen the fastening screw of the connector socket.
- ▶ Using pliers on the inner ring, pull the connector socket completely out of the ceiling rod with the cable.

- ▶ Before shortening the ceiling rod, remove the cable from the ceiling rod.
- ▶ Shorten the upper end of the ceiling rod to the desired length with a hacksaw and deburr it.

- ▶ Remove fixing screw "C".
- ▶ Insert the ceiling rod "A" into the ceiling bracket "B" and drill through the existing hole in the ceiling bracket with $d = 9\text{mm}$. Drill the opposite hole separately.
- ▶ Note: After sawing and drilling, pull the cable through from the bottom side of the rod to the top side of the rod (3-pin connector first)

- ▶ Insert the cable with the connector socket back into the ceiling rod.
- ▶ Align the threaded hole in the connector socket with the existing hole in the ceiling rod and fix it with the fastening screw.

3.3 Ceiling bracket installation

DANGER

Installation by qualified personnel

- ▶ Installation may only be carried out by a qualified specialist. Without the appropriate knowledge, life-threatening dangers may present themselves.
- ▶ Two people are required for the installation

DANGER

Life-threatening danger from falling light.


- ▶ The ceiling mount may only be attached to ceilings with a concrete strength class of B25 (C20/25) or greater.
- ▶ Reinforcement parts of the solid floor must not come into contact. In case of doubt, a licensed specialist must confirm the installation on the respective installation base. The load-bearing capacity of the ceiling construction must first be planned, checked and confirmed by a structural engineer.
- ▶ The drilling must be carried out professionally in compliance with the drilling tolerances approved by the manufacturer of the fastening anchor. In the event of incorrect drilling - for example when drilling into a rebar - a structural engineer must be called in.
- ▶ Install the light in such a way that the height stops are not permanently stressed during operation.
- ▶ In the case of plaster or cladding coating the concrete, the fastening anchor must be driven completely into the concrete.
- ▶ The screws must be carefully tightened using a torque spanner according to the manufacturer's instructions for the fixing anchors.

DANGER

Danger to life through an electric shock.


- ▶ The mains cable must be disconnected from the mains at all poles by an external lockable switch and secured against being switched on again.


- ▶ Draw in 6 drilling marks.
- ▶ Note the position of the opening "A" for the purpose of power connection

WARNING


Wear protective equipment according to the tool manufacturer's instructions

- ▶ Drill holes and blow out with bellows

- ▶ Check the distances between the drilled holes


- ▶ Hold the ceiling bracket to the ceiling and hammer in the fastening anchor.
- ▶ Tighten the fastener according to the manufacturer's instructions.

DANGER


Danger to life through an electric shock.

- ▶ Do not switch on the mains power supply before the light head has been installed
- ▶ To avoid the risk of an electric shock, this device may only be connected to a supply network with a protective earth conductor


- ▶ Create a network connection.


3.4 Ceiling bracket installation


- ▶ Pull the cable of the ceiling rod through the opening "B" of the ceiling bracket.
- ▶ Insert the ceiling rod "A" in the ceiling bracket.
- ▶ Secure with locking screw "C" and M8 nut.

- ▶ Tighten locking screw "C" M8 and nut (20 NM).
- ▶ Tighten all 4 setscrews "D" (5 NM).

- ▶ Connect the connector of the ceiling rod to the connector of the power unit.


- ▶ Slide the ceiling hood and ring over the ceiling bracket and screw tight.


3.5 Ceiling arm installation

WARNING


- ▶ **Risk of injury** from spring arm springing up
- ▶ The spring arm can open suddenly when the binding agent is removed and cause injuries. Please remove the binder carefully.


- ▶ Carefully turn and remove covers "A" from joint cover "B".
- ▶ Loosen the screws "C" and remove the joint covers "B" from the transom.


- ▶ Carefully remove the rotary connector "D".


- ▶ Insert the ceiling arm into the ceiling rod.

- ▶ After inserting the ceiling arm, first install the ring "A", then the locking ring "B" followed by the closing ring "C".

- ▶ Carefully insert the rotary connector "D".

- ▶ Install joint covers "B" and screw in screws "C".
- ▶ Insert covers "A" and secure by turning.

- ▶ To continue with the installation of the **light head**, see **Section 6**.

4. INSTALLATION: Triango 60 W

4.1 Load data

Bending moment M_B	275 Nm
Gravitational force F_G	135 N

4.2 Wall bracket installation

- **Fastening material** is **not included** in the scope of delivery.

DANGER

Installation by qualified personnel

- Installation may only be carried out by a qualified specialist. Without the appropriate knowledge, life-threatening dangers may present themselves.
- Two people are required for the installation

CAUTION

Use the load data table to determine the fasteners

- Note the dimensions of the linkage before installation

CAUTION


Note the position of the wall bracket

- The wall bracket must be aligned according to the axis in the picture
- Disregarding the correct alignment will result in mechanical safety defects
- For lightweight walls, we recommend a counter plate (not included in the scope of delivery)

DANGER

Danger to life through an electric shock.


- The mains cable must be disconnected from the mains at all poles by an external lockable switch.


- Draw in 6 drilling marks.
- Note the position of the opening "A" for the purpose of power connection.

WARNING


Wear protective equipment according to the tool manufacturer's instructions

- Drill holes and blow out with bellows.

- Check the distances between the drilled holes.


- Position the wall bracket on the wall and hammer in the fastening anchor.
- Tighten the fastener according to the manufacturer's instructions.


DANGER

Danger to life through an electric shock.


- Do not switch on the mains power supply before the light head has been installed.
- To avoid the risk of an electric shock, this device may only be connected to a supply network with a protective earth conductor

- Create a network connection.


- Insert wall bracket "A" (with end ring "C" and hood "B" attached) into wall bracket "D" and at the same time pull the connector through the rectangular recess.

- Install the wall bracket in vertical alignment with the locking screw and nut M8 "E" and tighten (20 Nm).
- Tighten all 4 setscrews "F" (5 NM).

- Connect the connector of the ceiling rod to the connector of the power supply unit.


- Slide the hood to the wall and screw tight with a ring (0.5Nm).

4.3 Wall arm installation


WARNING


- **Risk of injury** from spring arm springing up
- The spring arm can open suddenly when the binding agent is removed and cause injuries. Please remove the binder carefully.


- Pull off the covers "A" carefully.
- Unscrew screws "C" and remove joint covers "B".


- Carefully remove the rotary connector "D".


- Install the wall arm.

- After inserting the wall arm, first install the ring "A", then the locking ring "B" followed by the closing ring "C".

- Carefully insert the rotary connector "D".

- Install joint covers "B" and screw in screws "C".
- Attach covers "A".

DANGER


Danger to life through an electric shock.

- Do not switch on the mains power supply before the light head has been installed


- To continue with the installation of the light head, see **Section 6**.


5. INSTALLATION: Triango 60 F


- Remove the M3 screws at the sides and then remove the power supply unit tub "A".


- Lead the upright rod cable through the roller stand.
- Fasten the upright rod "B" to the roller stand using the 2 Allen screws and serrated washers "C" (10 Nm).
- The cable holder must be oriented **to the rear**.


- Pull off the covers "A" of the spring arm carefully.

- Remove the M4 hexagon socket screws together with the washer.

- ▶ Connect the connector of the spring to the connector of the upright rod.


- ▶ Place the spring arm onto the upright rod.


CAUTION

Do not push the connection cable back into the rod


- ▶ The connector can be damaged by squeezing the cable.


- ▶ If the connection cable, which runs through the joint, is too long, it can be pulled back through the upright rod into the roller stand.


- ▶ Pull the cable back out of the upright rod and roll it up in the roller stand.


- ▶ Connect the same coloured strands of the lower upright rod and the power supply unit.

- ▶ Mount the power supply unit tub again using M3 screws and serrated washers.

- ▶ Align the thread in the spring arm with the opening in the upright rod, then insert the M4 hexagon socket screw and washer.

- ▶ Attach the covers "A" again, one after the other. Ensure that the cover "B" is positioned in the groove of both covers "A".

DANGER

Danger to life through an electric shock.

- ▶ Do not switch on the mains power supply before the light head has been installed

6. INSTALLATION: Triango 60 light head

DANGER


Danger to life through an electric shock.

- ▶ When performing any work on the device, disconnect it from the power supply or pull the connector out of the socket and secure the device from being switched on again.


WARNING

Risk of injury

- ▶ The spring arm is under a high spring load. If the boom is lowered without the device installed, it **MUST** be held steady.
- ▶ If released, it will spring up, which can cause serious injury.
- ▶ Only remove the body of the light when the boom is in the up position or when a second person is securely holding it in a lowered position.
- ▶ Always have a second person to help when installing or removing the body of the light, to avoid serious injury or damage.

- ▶ Connect the connector of the light head to the connector of the spring arm.
- ▶ Push the light head onto the spring arm.
- ▶ When aligning the light head, make sure that the flat side of the brass part points down to "A".

- ▶ Screw the grub screw into the locking part "B"
- ▶ The screw grips the previously aligned, flat side of the brass part

WARNING

Risk of injury from falling light head.

- ▶ Make sure the safety element is properly installed.
- ▶ Material damage due to improperly installed light head.
- ▶ Tighten the screw as required by the friction of the light head

- ▶ Attach handle.

7. OPERATION

Triango 60 C, W, F

WARNING

- ▶ This product emits possibly hazardous optical radiation. Do not stare at the light emitted from the surgical luminaire. Eye injury may occur.
- ▶ The optical radiation emitted by this product complies with the exposure limits for reducing the risk of photobiological hazards in IEC 60601-2-41.

Triango 60-1 C, W, F

- ▶ Switch the light on or off with button "1"
- ▶ Continuously dim the light by pressing and holding down button "2".

Triango 60-3 C, W, F

- ▶ Switch the light on or off with button "3"
- ▶ Continuously dim the light by pressing and holding down button "3".
- ▶ The light colour can be set using button "4" (3700K / 4300K / 4700K)

Triango Fokus 60-1 / 60-3 C, W, F


- ▶ Switch the light on or off with button "3"
- ▶ Continuously dim the light by pressing and holding down button "3".
- ▶ The focus can be set using button "5" (18 cm / 23 cm / 25 cm)

Triango 60 F**⚠ DANGER****Danger to life through an electric shock.**

- ▶ Do not plug in damaged power cables
- ▶ If there is any sign of damage to the power cable, replace it with a new one immediately
- ▶ Connection voltage and frequency must match the data on the type plate.
- ▶ Only connect to the supply network with a protective conductor

CAUTION**When changing location**


- ▶ Position the end device in the lowest position

- ▶ Unlock casters
- ▶ Do not run over any objects or the connection cable
- ▶ Hold tight during transport within the hospital.
- ▶ Beware of inclines, door sills, unevenness or other obstacles

CAUTION

- ▶ When not in use, wrap the power cable in the cable holders

- ▶ Plug in the power cable
- ▶ Connect the cable to the mains

8. CLEANING AND DISINFECTION**⚠ DANGER****Danger to life through an electric shock**

- ▶ Before performing disinfection cleaning, pull the connector out of the socket and secure the device against unintended switching on again. (new)

CAUTION**Damage caused by improper cleaning**

- ▶ Only such agents may be used for cleaning that do not impair the functionality of the light.
- ▶ No cleaning agents containing solvents, chlorine or abrasives may be used for cleaning, as such agents can lead to cracking in plastic parts, among other things.
- ▶ The agents used must be approved for use with plastics such as PC, PMMA, PA and ABS.
- ▶ Damage to the light due to concentrated disinfectants.
- ▶ For concentration and exposure time, please note the information on the supplementary sheet for the agent used.
- ▶ Scratches from unsuitable cloths.

RECOMMENDED DISINFECTANTS

- ▶ Lysoformin ®
- ▶ Dismozon ®
- ▶ Hexaquart ®plus
- ▶ Sagrotan ® rapid disinfectant cleaner

CAUTION**Dirt reduces the luminosity**

- ▶ Keep the screen clean by cleaning it regularly
- ▶ Wipe clean only

- Clean the PA screen with a leather cloth soaked in window cleaner.

CAUTION

In order to minimise the risk of disease transmission, the applicable occupational health and safety regulations and the requirements of the national responsible committees for hygiene and disinfection must be observed in addition to these operating instructions.

8.1 Sterilisation of the handle

- During sterilisation, **ISO 17665-1** (sterilisation of medical devices in moist heat) must be observed.

CAUTION

Damage to the handle

- Do not sterilise with hot air
- The handle must be packed in a sterile bag before sterilisation
- The handle is designed exclusively for steam sterilisation with triple fractionated pre-vacuum and saturated steam with the following parameters:

Temperature	134°C
Excess pressure	2.0 bar
Retention time	6 min
Drying in vacuum	20 min

- Check the mechanical integrity of the handle after each sterilisation
- Stop using the handle if it is damaged

9. SAFETY CONTROLS

DANGER

Danger to life through an electric shock.

- Disconnect the plug from the mains.
- The connecting cable must be checked for damage at least once a year.

CAUTION

- Maintenance and repairs may only be carried out by qualified electricians.
- The corresponding user profile is in Section 2 Safety Instructions.

ANNUALLY:

- Check the connecting cable for damage and replace if necessary
- Check for paint damage/cracks on plastic parts
- Check for deformation or damage to the support system
- Check for part loosening

9.1 SETTING SPRING FORCE


CAUTION

- The spring force is ideally set at the factory

Triango 60 C / Triango 60 W

CAUTION

- The lamp head **must** be installed before the spring force is adjusted.


- Move the spring arm into the highest possible position.
 - Insert a 6 mm Allen key through the hole in the central joint and adjust the spring force by turning the screw.
 - Turn the screw clockwise (+) to increase the spring force (as the spring arm lowers).
- Turn the screw anticlockwise (-) to decrease the spring force (as the spring arm rises).


Triango 60 F

CAUTION

- The lamp head **must** be installed before the spring force is adjusted.

- Move the spring arm into the highest position.
- Remove the side plastic covers and push the round plastic flap up along the spring arm

- Adjust the screw with a 4mm Allen key
- Turn the screw clockwise (+) to increase the spring force (as the spring arm lowers)
Turn the screw anticlockwise (-) to decrease the spring force (as the spring arm rises)

10. DEINSTALLATION

DANGER

Danger to life through an electric shock.

Before disassembly, the light must be completely disconnected from the mains voltage.

WARNING

Risk of injury

The spring arm is under a high spring force. If the end device is not removed in the uppermost spring arm position, the spring arm will spring up and can cause serious injuries. Only disassemble the end device when the spring arm is in the uppermost position

10.1 Disposal

Do not put the light in the household waste. Dispose of the light at a disposal point in accordance with local regulations or give it to a dealer with the appropriate service.
Cut the cable directly at the housing.

The products listed above are over 95% recyclable. The lights are designed to be recyclable so that a high proportion of the materials used can be recycled or used to generate energy after the end of the service life of these products. They do not contain any dangerous substances or substances that require monitoring.

11. ACCESSORIES

- Handle (item number D10.442.000)

- Grip cover (item number D15.445.000)

12. ADDITIONAL INFORMATION

The light itself is maintenance-free.

Additional documentation for this product can be requested from the manufacturer on request.

The use of this light does not create any risks that could affect other devices.

To save energy, the light should only be switched on when it is actually being used.

All serious incidents related to the device **must be reported** to the manufacturer or its representative and to the competent authority of the Member State where the user is established.

13. TROUBLESHOOTING

Fault	Possible cause	Troubleshooting	User profile
Light does not come on	Contact fault	Switch on again	All
Light does not come on	Light defective	Contact manufacturer	Only via manufacturer
Light does not come on	No mains voltage	Check voltage, inspect all connections	Electrician

14. TECHNICAL DATA

Electrical values:	
Rated connection voltage	100-240 V
Frequency range	50-60 Hz
Power consumption:	
Triango 60-1 / Triango 60-3	21 W (max. 45VA)
Triango Fokus 60-1	33 W (max. 64VA)
Integrated electronic transformer	24VDC output
Photometric values*:	
Central illuminance Ec at 1.0m distance	60,000 lx
Illuminated field diameter d10 at a distance of 1.0m:	
Triango 60-1 / Triango 60-3	Ø = 18 cm
Triango Fokus 60-1	Ø = 18 cm / 23 cm / 25 cm
Illuminated field diameter d50 at a distance of 1.0m	Ø = 10 cm
Colour temperature:	
Triango 60-1 / Triango Fokus 60-1	4700K
Triango 60-3	3700K / 4300K / 4700K
Colour rendering index Ra	95
Residual illuminance when shaded by a shadow	< 1%
Residual illuminance when shaded by two shadows	64.9%
Residual illuminance in tube	100%
Residual illuminance in tube with a shadow	< 1%
Residual illuminance in tube with two shadows	65.1%
Illumination depth L1 + L2	116 cm
* -10% / +20% tolerance	

*	TRIANGO 60-1	TRIANGO 60-3	TRIANGO Fokus 60-1
Maximum illuminance distance (D _{MI})	200 mm	200 mm	200 mm
Maximum irradiance and radiance settings	200 mm max. illuminance	200 mm @4700K max. illuminance	200 mm max. illuminance wide focus position 25 cm
Total irradiance	268,7 W/m ²	456,8 W/m ²	551,3 W/m ²
Light field diameter d ₁₀	Measured d ₁₀ : 184 mm	Measured d ₁₀ : 182 mm	Measured d ₁₀ : 186 mm 226 mm 246 mm
ratio d ₅₀ /d ₁₀	0.57	0.58	0.55 0.50 0.50
Colour rendering index R9	76.4	90.2	80.7
* -10% / +20% tolerance			

Ambient conditions for transport, storage and operation:

Ambient temperature (storage and transport)	-20°C to +70°C
Ambient temperature (operation)	+10°C to +35°C
rel. Humidity (non-condensing) (operating)	max. 75%
Maximum use height (operating) standard version (H) version	3000 m (above sea level) 5000 m (above sea level)
Atmospheric pressure standard version (H) version	70-106 kPa 54-106 kPa

Weight:

Light head	1.5 kg
Triango 60 C	12.5 kg
Triango 60 W	12.0 kg
Triango 60 F	17 kg

Type of operation

Type of operation	Continuous operation
-------------------	----------------------

Classification

Triango 60 C / W / F	Protection class I
Type of protection	IP 20
Light head	IP 43 (horizontal position)
Classification according to EU REGULATION 2017/745 (MDR), Item 51 U.S. FDA Device Class	Class I Class I
Electrical safety test and EMC according to:	AAMI ES60601-1: 2005/A2:2021 CAN/CSA-C22.2 No 60601-1:14 IEC 60601-1:2005 + A2:2020 IEC 60601-1-2:2014 + A2:2020 IEC 60601-2-41:2021
Blue light hazard according to IEC 62471:2006;modified	RG2 (medium risk)

Lifetime of the light source:

Lifetime

50,000h (L70/B50)

15. NAMING THEME

TRIANGO (X) 60-(Y) (Z) (W) (V) (U)

(U): **CH, USA, SCHUKO, AUS, UK, BR** or **US-HG**. country versions(V): **H** or not. H means higher max operating use height of 5000m(W): **T** or not. T means use of the Derungs Teno arm system.(Z): **C** or **W** or **F**. C= ceiling version, W= wall version, F= mobile stand version(Y): **1** or **3**. 1= light single-colour version, 3= light three-colour version

(60): Product Family 60, stands for Ec 60Klux @1m distance

(X): **Fokus** or not, Focus stands for the light adjustment function

Product Family Name

16. ELECTROMAGNETIC COMPATIBILITY (EMC)

Medical electrical devices are subject to special precautionary measures with regard to electromagnetic compatibility. This device can be affected by other electrical devices.

WARNING

Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the Triango 60, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

Electromagnetic environment

The device may only be operated in the environments specified in the "Intended Use" section of the Operating Instructions. The medical device is intended for use in the electromagnetic environment specified below.

The above listed models are intended for use in the electromagnetic environment specified below. The customer or the user of the above listed models should assure that they are used in such an environment.		
Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The above listed models use RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations / flicker emissions IEC 61000-3-3	Complies	

The above listed models are intended for use in the electromagnetic environment specified below. The customer or the user of the above listed models should assure that they are used in such an environment.			
Enclosure Port			
Immunity test	Test Condition	IEC 60601 Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV Contact $\pm 2,4,8,15$ kV Air	± 8 kV Contact $\pm 2,4,8,15$ kV Air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Radiated RF EM fields and Proximity fields from RF wireless communications equipment IEC 61000-4-3	10 V/m 80 MHz – 2,7 GHz 80% AM 1kHz	10 V/m 80 MHz – 2,7 GHz	Mains power quality should be that of a professional healthcare facility environment and Home healthcare environment.
	385MHz (18Hz Pulse Modulation)	27 V/m	
	450MHz (FM+/-5KHz deviation 1kHz sine or 18Hz Pulse Modulation)	28 V/m	
	710MHz (217Hz PM)	9 V/m	
	745MHz (217Hz PM)	9 V/m	
	780MHz (217Hz PM)	9 V/m	
	810MHz (18Hz PM)	28 V/m	
	870MHz (18Hz PM)	28 V/m	
	930MHz (18Hz PM)	28 V/m	
	1720MHz (217Hz PM)	28 V/m	
	1845MHz (217Hz PM)	28 V/m	
	1970MHz (217Hz PM)	28 V/m	
	2450MHz (217Hz PM)	28 V/m	
	5240MHz (217Hz PM)	9 V/m	
	5500MHz (217Hz PM)	9 V/m	
	5785MHz (217Hz PM)	9 V/m	
RATED power frequency magnetic fields IEC 61000-4-8	50Hz or 60Hz	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
IMMUNITY to proximity magnetic fields	30 kHz CW Test level: 8 A/m 134,2 kHz PM 2,1 kHz Test level: 65 A/m 13,56 MHz PM 50 kHz Test level: 7,5 A/m	Not applicable.	The Triango 60 light do not contain magnetically sensitive components or circuitry within the ENCLOSURE.

The above listed models are intended for use in the electromagnetic environment specified below. The customer or the user of the above listed models should assure that they are used in such an environment.

Input a.c. power PORT

Immunity test	Test Condition	IEC 60601 Compliance level	Electromagnetic environment - guidance
Electrical fast transient/bursts IEC 610004-4	± 2 kV 100kHz Repetition frequency	± 2 kV 100kHz Repetition frequency	Mains power quality should be that of a professional healthcare facility environment and Home healthcare environment.
Surges IEC 61000-4-5	$\pm 0,5$ kV, ± 1 kV line(s) to line(s) $\pm 0,5$ kV, ± 1 kV , ± 2 kV line(s) to ground	± 1 kV , Differential mode ± 2 kV Common mode	Mains power quality should be that of a professional healthcare facility environment and Home healthcare environment.
Conducted RF induced by RF fields IEC 61000-4-6	3 V 0,15 MHz – 80 MHz 6 V m) in ISM and amateur radio bands between 0,15 MHz and 80 MHz (n) 80 % AM at 1 kHz	3 V 0,15 MHz – 80 MHz 6 V m) in ISM and amateur radio bands between 0,15 MHz and 80 MHz (n) 80 % AM at 1 kHz	Mains power quality should be that of a professional healthcare facility environment and Home healthcare environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0% U_T ; 0°,45°,90°,135°,180°,225°,270°,315° 0% U_T ; 0° 0% U_T ; 70% 0% U_T ; 0%	0,5 Cycles 1Cycle 25/30 Cycles (50/60Hz) 250/300 Cycles (50/60Hz) (5s)	Mains power quality should be that of a professional healthcare facility environment and Home healthcare environment. If the user of the above listed models requires continued operation during power mains interruptions, it is recommended that the above listed models are powered from an uninterruptible power supply or battery.

Comment: n) The ISM (industrial, scientific and medical) bands between 0,15 MHz and 80 MHz are 6,765 MHz to 6,795 MHz; 13,553 MHz to 13,567 MHz; 26,957 MHz to 27,283 MHz; and 40,66 MHz to 40,70 MHz. The amateur radio bands between 0,15 MHz and 80 MHz are 1,8 MHz to 2,0 MHz, 3,5 MHz to 4,0 MHz, 5,3 MHz to 5,4 MHz, 7 MHz to 7,3 MHz, 10,1 MHz to 10,15 MHz, 14 MHz to 14,2 MHz, 18,07 MHz to 18,17 MHz, 21,0 MHz to 21,4 MHz, 24,89 MHz to 24,99 MHz, 28,0 MHz to 29,7 MHz and 50,0 MHz to 54,0 MHz.

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USA

Waldmann Lighting Company
9, W. Century Drive
Wheeling, Illinois 60090
USA
Telephone +1 847 520 1060
Fax +1 847 520 1730
www.waldmannlighting.com
waldmann@waldmannlighting.com

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Hofmattstrasse 12
CH-9200 Gossau SG
Switzerland
Telephone +41 71 388 11 66
Fax +41 71 388 11 77

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

LS, Produktmanagement, 08.10.2024


Installation Height

Installation Height (with intermediate ceiling adapter)

Installation Position

* If installation height = 195cm,
and bed height = 80cm
and lighting distance = 100cm (80'000 lux)

**If possible, we recommend that the light is
installed opposite the treating doctor.
This makes operating the light more ergonomic.

Ceiling base plate – Hole pattern

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Derungs Medical GmbH