Montage- und Gebrauchsanweisung	04	DEU
Installation and operating instructions	22	ENG
Mode d'emploi	40	FRA
Istruzioni per il montaggio e l'uso	58	ITA
Instrucciones de montaje y funcionamiento	76	ESP
Instruções de instalação e de utilização	94	POR
Montage- en gebruiksaanwijzing	112	NLD
Bruksanvisning	130	SWE
取り付けおよび使用の説明書	148	JAP

D ^{med} ®	Triango	100 C
D ^{med} ®	Triango	100 W
D ^{med} ®	Triango	100 F

Operationsleuchte
 Surgical luminaire
 Luminaire chirurgical
 Lampada chirurgica
 Luminaria quirúrgica
 Candeeiro de cirurgia
 Operatiearmatuur
 Kirurgisk armatur
 手術用照明器具

SYMBOLE; SYMBOLS; SYMBOLES; SIMBOLI; SÍMBOLOS; SÍMBOLOS; SYMBOLEN; SYMBOLER; 記号の意味

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
IMPORTANT!
THIS OPERATING MANUAL MUST BE THOROUGHLY READ
BEFORE USING THE PRODUCT!
→ STORE FOR LATER USE!

CONTENT

1.2	TRIANGO 100 W	23
1.3	TRIANGO 100 F	23
2.	SAFETY INSTRUCTIONS	24
2.1	Intended use	24
2.2	User profiles	24
2.3	Safety instructions	24
2.4	Warning levels	24
2.5	Specific fastening instructions	25
3.	INSTALLATION: TRIANGO 100 C	25
3.1	Workload data	25
3.2	Shortening and drilling the ceiling tube	25
3.3	Installing the ceiling holder	25
3.4	Installing the ceiling tube	26
3.5	Installing the ceiling arm	27
4.	INSTALLATION: TRIANGO 100 W	27
4.1	Workload data	27
4.2	Installing the wall holder	27
4.3	Installing the wall mounting arm	28
5.	INSTALLATION: TRIANGO 100 F	29
6.	INSTALLATION: luminaire head	30
7.	OPERATION	31
7.1	Operating function	31
7.2	Explanation of light modus.....	31
8.	CLEANING AND DISINFECTION	32
8.1	Sterilizing the hand grip	33
9.	SAFETY INSPECTIONS	33
9.1	Setting the spring force.....	33
10.	DISMANTLING	33
10.1	Disposal.....	34
11.	ACCESSORIES	34
12.	ADDITIONAL INFORMATION	34
13.	TROUBLESHOOTING	34
14.	TECHNICAL DATA	35
15.	ELECTROMAGNETIC COMPATIBILITY (EMC)	37

1. VERSIONS AND DELIVERY CONTENT


1.1 TRIANGO 100 C

A: Ceiling bracket, ceiling bonnet and end ring
C: Lamp head with sterilisable handle

B: Ceiling tube
D: Ceiling arm


1.2 TRIANGO 100 W

A: Wall bracket, wall bonnet and end ring
C: Wall arm

B: Wall joint
D: Lamp head with sterilisable handle

1.3 TRIANGO 100 F

A: Lamp head with sterilisable handle
C: Lower standtube

B: Upper stand tube with spring arm
D: Rolling frame
E: Mains cable

2. SAFETY INSTRUCTIONS

2.1 Intended use

The TRIANGO 100 light is a surgical luminaire. It is a single light in the patient environment for use in treatment rooms to support diagnosis or treatment that does not pose a hazard to the patient in the event of interruption due to a light failure. It is intended for continuous operation and is not intended to be combined with other medical devices.

The essential performance of TRIANGO 100 is the delivery of sufficient light to the operating field while limiting the delivery of radiant energy to the operating field and its observer. This is achieved by providing the capability to deliver central illuminance of at least 40 Klx, while ensuring the total irradiance at the light field centre at the maximum illuminance distance does not exceed 700 W/m².

2.2 User profiles

Health Professional

All individuals who have completed medical training and work in the professional field they trained for.

Cleaning specialist

Trained in national and workplace hygiene regulations.

Qualified electrician

Trained in electronics and electrical technology and knows the relevant standards and regulations.

Qualified specialist

Qualified due to his technical training, knowledge and experience and knowledge of the rules, to perform assembly/dismantling.

2.3 Safety instructions

- ▶ Operation by health professionals
- ▶ This manual is part of the product and must be stored and made available to all future users.
- ▶ All work on the luminaire (including repairs) may only be performed by a qualified electrician. Assembly may only be done by a qualified specialist.
- ▶ The luminaire may not be altered or tampered with. Only approved original parts may be used. Any use other than that intended, using original parts, can change technical parameters and pose a death hazard.
- ▶ Do not exceed the maximum weight and do not hang, lean or climb on the arm, because this may cause the device to tip over, which could result in serious injury.
- ▶ Operation in potentially explosive areas is prohibited. The luminaire's current supply is a potential ignition source.
- ▶ The luminaire must be used only in dry, dust-free rooms.
- ▶ The luminaire should not be left on without supervision.
- ▶ Only connect the luminaire to the supply network with a protective ground conductor (PE), to prevent electrical shock.
- ▶ In the case of luminaires of protection class I, the protective ground conductor (PE) must be connected to luminaires housing.
- ▶ Do not use a damaged luminaire. Defective cables and a defective hand grip also pose a potential hazard. Do not place a cord near heat sources or sharp edges.
- ▶ Do not place extra loads on the luminaire head and arm system.

- ▶ The luminaire must not be covered with a cloth or similar item during operation.
- ▶ The ventilation openings (if present) must always be kept clear during operation!
- ▶ The luminaire must not be operated near external heat sources that exceed the luminaire's maximum ambient temperature.
- ▶ The luminaire must not be used outside the specified ambient conditions.
- ▶ Do not use with medical devices that may react sensitively to a light spectrum within the visible range (such as pulsating light and/or light with high illumination intensity).
- ▶ This product emits possibly hazardous optical radiation. Do not stare at the light emitted from the luminaire. Eye injury may occur.
- ▶ The luminaire may only be used for the intended use described here.
- ▶ The manufacturer cannot be held responsible for any damages resulting from use deviating from the intended use, or failure to observe the safety instructions and warnings.
- ▶ When using more than one luminaire at the same time, the total illumination intensity must not exceed E_e 1000W/m².
- ▶ Before connecting to the power grid, it is essential to check that the grid data match the device data.
- ▶ For operations, only the spot function with an illumination field of 18 cm should be used with the Fokus version. Further focus options (illumination field 23 cm or 28 cm) are only additional functions for support.
- ▶ The light must be subjected to an annual safety inspection. This is described in chapter 9.
- ▶ **TRIANGO 100 F**
The light must be held securely during transport within the hospital.

2.4 Warning levels

DANGER

Warning of hazards that can result in **death or serious injury** if instructions are not followed.

WARNING

Warning of hazards that can result in **injury** if instructions are not followed.

CAUTION

Warning of hazards that can cause **material damage** if instructions are not followed.

2.5 Specific fastening instructions

TRIANGO 100 C

- ▶ **Fastening material is not included with delivery.**
- ▶ The ceiling must only be mounted on ceilings having a concrete stability class B25 (C20/25) or higher.
- ▶ Solid ceiling reinforcement parts must not come into contact with the ceiling mounting. If there is any doubt, an authorized professional must confirm that mounting is possible on the mounting base being used. The load-bearing capacity of the ceiling construction must be planned, checked and confirmed by a structural engineer beforehand.
- ▶ Borings must be made by specialists in compliance with the boring tolerances for the reinforcement anchor that have been approved by the manufacturer. In the event of faulty drilling, such as drilling through a metal support bar, a structural engineer must be called in.
- ▶ The luminaire should be installed in such a way as to avoid permanent strain on the height stops during operation.
- ▶ If there is plastering or cladding over the concrete, make sure the fixing anchor is properly attached right into the concrete.
- ▶ Screws must be tightened carefully using a torque wrench that complies with the instructions of the mounting element manufacturer.

TRIANGO 100 W


- ▶ **Fastening material is not included with delivery.**
- ▶ The luminaire must only be mounted on walls that can guarantee a secure hold. Specialist staff will find details of the requirements in Chapter 5.1 (Load data).

3. INSTALLATION: TRIANGO 100 C


3.1 Workload data

Bending moment M_B	245 Nm
Vertical weight F_G	160 N

3.2 Shortening and drilling the ceiling tube

- ▶ Before shortening the ceiling tube, remove the cable from the ceiling tube.
- ▶ Use a hack saw to shorten the upper end of the ceiling tube to the desired length, and then deburr it.

- ▶ Remove the fastening screw "C".
- ▶ Insert the ceiling tube "A" into the ceiling bracket "B" and drill the ceiling holder's existing hole to 9 mm in diameter. Drill the opposite hole separately.
- ▶ Note: Pull the cable through from the lower side to the upper tube side (the 3-pin plug first) after sawing and drilling.

3.3 Installing the ceiling holder

DANGER

Assembly by qualified personnel

- ▶ Assembly must be done by qualified personnel only. Lack of appropriate knowledge may pose a death hazard.
- ▶ Two people are needed for assembly.

DANGER


Death hazard from falling luminaire.

- ▶ The ceiling mount must only be mounted on ceilings with a concrete stability class of B25 (C20/25) or higher.
- ▶ There should be no contact with reinforcement parts of the solid ceiling during mounting. When in doubt, a licensed technician must approve installation on the specific installation surface. The load bearing capacity of the ceiling structure must be planned, checked and confirmed by a structural engineer beforehand.
- ▶ The holes must be drilled by a professional and in compliance with the drilling tolerances permitted by the manufacturer of the fastening anchor. If a drilling error occurs, such as drilling into a reinforcement bar, a structural engineer must be consulted.
- ▶ Install the luminaire such that the vertical stops are not continuously stressed during operation.
- ▶ With plaster or cladding in front of concrete, the fastening anchor must be driven completely into the concrete.
- ▶ The screws must be tightened carefully using a torque wrench according to the specifications of the fastening anchor manufacturer.

DANGER

Danger to life through an electric shock.


- ▶ The mains cable must be disconnected from the mains at all poles by an external lockable switch and secured against being switched on again.


- Draw 6 drill marks.
- Observe the position of opening "A" for connecting the current.

WARNING


Wear safety equipment according to the tool manufacturer's instructions.

- Drill holes and blow out with a bellows.

- Check the hole distances.


- Hold the ceiling bracket to the ceiling and hammer in the fastening anchor.
- Tighten the fastener according to manufacturer's instructions.

DANGER

Death hazard from electric shock.

- Do not turn on the power supply before the luminaire head has been installed.
- This device must only be connected to a power grid with a protective ground conductor (PE) to prevent risk of electric shock

- Connect the power.


3.4 Installing the ceiling tube

DANGER


Risk of death from falling luminaire.

- The following safety screw "C" and the M8 nut, as well as the grub screws "D" are safety-relevant components and must always be installed.


Before installing the ceiling tube, the two holes in the ceiling tube for the safety screw "C" must always be drilled (see chapter "3.2 Shortening and drilling the ceiling tube")


- Pull the ceiling tube cable through the opening "B" of the ceiling holder.
- Insert the ceiling tube "A" into the ceiling holder.
- Secure the safety screw "C" and the M8 nut.

- Tighten the M8 safety screw "C" and the nut (20 Nm).
- Tighten all 4 grub screws "D" (5 Nm).

- Connect the ceiling tube plug with the power supply plug.


- Slide the ceiling cover and ring over the ceiling bracket and screw them tight.

3.5 Installing the ceiling arm

WARNING


The opening spring balanced arm is an injury hazard

- ▶ Risk of injury from spring arm springing up
- ▶ The spring arm can open suddenly when the binding agent is removed and cause injuries. Please remove the binder carefully.


- ▶ Connect the plug from the ceiling tube with plug on the ceiling arm
- ▶ Carefully push protruding cable into the ceiling tube

Perform this assembly step on an even surface:

- ▶ Insert the sliding piece 1 above into the ceiling arm groove 2
- ▶ Push the ceiling tube 3 and ceiling arm 2 together
- ▶ Secure the sliding piece 1 by screwing in the threaded stud 4 (up to the stop)
- ▶ Loosen the threaded stud 4 again by a maximum of a ¼ turn and checking the rotating function.

- ▶ Continue with luminaire head assembly, see Chap. 6

4. INSTALLATION: TRIANGO 100 W

4.1 Workload data

Bending moment M_B	275 Nm
Vertical weight F_G	155 N

4.2 Installing the wall holder

- ▶ **Fastening hardware** does **not** come with delivery.

DANGER

Assembly by qualified personnel

- ▶ Assembly must be done by qualified personnel only. Lack of appropriate knowledge may pose a death hazard.
- ▶ Two people are needed for assembly.

CAUTION

Determine the fasteners to use from the workload data table

- ▶ Before assembly, make sure the rod dimensions are correct.

CAUTION


Make sure the wall holder is in the correct position.

- ▶ The wall holder must be aligned according to the axis in the image.
- ▶ Misalignment will compromise mechanical safety.
- ▶ We recommend using a counterplate on lightweight construction walls (not included with delivery).

DANGER

Death hazard from electric shock.


- ▶ All poles of the luminaire must be disconnected from the power supply by an external lockable switch.


- ▶ Draw 6 drill marks.
- ▶ Observe the position of opening "A" for connecting the current.

⚠ WARNING


Wear safety equipment according to the tool manufacturer's instructions.

- Drill holes and blow out with a bellows.

- Check the hole distances.


- Position the wall bracket on the wall and hammer in the fastening anchor.
- Tighten the fastener according to manufacturer's instructions.


⚠ DANGER

Death hazard from electric shock.


- Do not turn on the power supply before the luminaire head has been installed.
- This device must only be connected to a power grid with a protective ground conductor (PE) to prevent risk of electric shock

- Connect the power.


- Insert wall bracket "A" (with end ring "C" and cover "B" in position) into the wall holder "D" and simultaneously pull the plug through the rectangular cutout.

- When the wall bracket is vertically aligned, install it with the M8 locking screw and nut "E" and tighten (20 Nm).
- Tighten all 4 grub screws "F" (5 Nm).

- Connect the ceiling tube plug with the power supply plug.


- Push the cover against the wall and screw tight with the ring (0.5 Nm)

4.3 Installing the wall mounting arm


⚠ WARNING

The opening spring balanced arm is an injury hazard


- Risk of injury from spring arm springing up.
- The spring arm can open suddenly when the binding agent is removed and cause injuries. Please remove the binder carefully.


- Connect the wall-mounted joint plug to the wall bracket plug


- Connect the wall arm plug to the wall-mounted joint plug

- Insert the arm for wall-mounting **3** into wall bracket **1** slightly
- Insert the sliding piece **2** into the groove
- Push up the arm **3** in such a way that the hole in the wall bracket is aligned with the hole in the sliding piece
- Screw grub screw **4** into wall bracket


- In the case of a dead stop, screw this back by a ¼ turn and check the turning function


- Continue with luminaire head assembly, see Chap.6

5. INSTALLATION: TRIANGO 100 F


- Remove the M3 screws on the sides and remove the power supply tray "A".


- Loosen the 2 screws on the lower support tube.
- Feed the support tube cable through the roller stand.
- Attach the support tube to the roller stand using the 2 Allen screws and serrated washers (10 Nm).
- The cable holder must be aligned rearward.


- Connect the same-colored stranded wires of the lower support tube and of the power supply.


- Reinstall the power supply tray with the M3 Allen screws and serrated washers.


- Loosen the grub screw **1**
- Place the spring balanced arm **2** on the upper support tube
- Attention: Grub screw **1** must be inserted into slot **3**.

- Screw in the grub screw in the alignment shown
- In the case of a grub screw dead stop, screw this back by a ¼ turn and check the turning function (Swing range $\pm 30^\circ$)

- Push together the plugs on the upper and lower support tube until they click.
- Place the upper support tube on the lower one.

- Tighten the cable holder using a type 3 Allen screw on the upper support tube (2.4 Nm).

WARNING

Never unscrew the upper cable holder — injury hazard.

- If both cable holders are unscrewed, the connecting piece will come loose and fall down, which can cause personal injury and damage to the cable and device.

DANGER

Death hazard from electric shock.

- Do not turn on the power supply before the luminaire head has been installed.

6. INSTALLATION: luminaire head

DANGER


Death hazard from electric shock.

- Before any work on the device, disconnect the current supply or pull the plug out of the socket and secure the device against being turned back on.


WARNING

Injury hazard


- The spring balanced arm is under high spring stress. If the boom is lowered without the installed device, it **MUST** be held firmly.
- If it is let go, it springs upward and can cause severe injury.
- Never remove the luminaire body unless the boom is in the upper position or is held securely in low position by a second person.
- To prevent severe injury or damage, always get the help of a second person when installing and removing the luminaire body.

- Connect the luminaire head plug to the wall-mounted joint plug

- Before assembly one must ensure that the brass nose (see arrow) and the cardan joint bracket are on the left


- Push together the luminaire head **1** and arm **3**; in doing so the cardan joint bracket must be on the left **2** of the luminaire head
- then screw in both countersunk screws M4 (both with a toothed washer) in (1.5 Nm)
- Check the turning function

⚠ WARNING

Injury hazard from falling luminaire head.

- ▶ Make sure the locking device is properly installed.
- ▶ Damage to due to a poorly mounted luminaire head.
- ▶ Tighten the screw as the luminaire head's friction requires.

- ▶ Attach the hand grip.


7. OPERATION

CAUTION

- ▶ This product emits potentially hazardous optical radiation. Do not stare directly into the light cone. Eye irritation may occur.
- ▶ The optical radiation emitted by this product complies with the exposure limits for reducing the risk of photobiological hazards in IEC 60601-2-41.


7.1 Operating function

The different functions depend on the model


- ▶ Switch the luminaire on or off using the button 1
- ▶ The luminaire can be dimmed using the +/- buttons 2
- ▶ The light color can be adjusted using button 3
- ▶ The focus can be adjusted using button 4
- ▶ The endoscope mode can be adjusted using button 5
- ▶ Before each use, perform a function test: all LEDs in the light cone must come on.


7.2 Explanation of light modus


- ▶ only LED 1 lights up: 40% illuminance
- ▶ LED 1 and 2 light up: 60% illuminance
- ▶ LED 1, 2 and 3 light up: 80% illuminance
- ▶ LED 1, 2, 3 and 4 light up: 100% illuminance

- ▶ only LED 1 lights up (blue): colour temperature = 3700 K
- ▶ only LED 2 lights up (white): colour temperature = 4300 K
- ▶ only LED 3 lights up (red): colour temperature = 4700 K

- ▶ only LED 1 lights up: illuminancefield @ 100 cm = 28 cm
- ▶ only LED 2 lights up: illuminancefield @ 100 cm = 23 cm
- ▶ only LED 3 lights up: illuminancefield @ 100 cm = 18 cm

- ▶ LED 1 lights up: ENDO mode is activated: 10% of the specified illuminance (according to chapter 14)

TRIANGO 100 F


⚠ DANGER

Electric shocks are life-threatening.

- ▶ Do not plug in any damaged power cables.
- ▶ If there are any signs of damage to the power cable, immediately replace it with a new one.
- ▶ Connection voltage and frequency must match data on the type plate.
- ▶ Only connect the luminaire to the supply network with a protective earth conductor (PE).

⚠ DANGER**In the event of a change of location**


- ▶ Put the terminal into its lowest position.

- ▶ Unlock the castors.
- ▶ Take care not to run over anything, including the connection cable.
- ▶ Hold firmly during transport within the hospital.
- ▶ Take particular care on slopes, thresholds, bumps and other obstacles.

CAUTION**TRIANGO 100 F**

- ▶ Wind the power cable onto the cable holder when not in operation

- ▶ Plug in cable
- ▶ Connect cable to the network
- ▶ Lock the castors

8. CLEANING AND DISINFECTION**⚠ DANGER****Death hazard from electric shock**

- ▶ Before disinfection cleaning, pull the plug out of the socket and secure the device against being turned back on.

CAUTION**Material damage due to incorrect cleaning**


- ▶ For cleaning, use only products that do not impair operation of the luminaire.
- ▶ For cleaning, do not use solvent-based, chlorine-based or abrasive detergents, because they can crack plastic parts and cause other damage.
- ▶ The cleaning agents must be approved for use on plastics such as PC, PMMA, PA and ABS.
- ▶ Concentrated disinfectant can damage the luminaire.
- ▶ For concentration and application times, check the information provided with the product used.
- ▶ The wrong cloth may cause scratches.

RECOMMENDED DISINFECTANTS


- ▶ Bacillol AF
- ▶ Cleanisept Wipes (Maxi)
- ▶ Descogen Oxy Wipes
- ▶ Descosept Sensitive Wipes
- ▶ Incidin Alcohol Wipe
- ▶ Incidin OxyWipe NG
- ▶ Incidin Plus
- ▶ Meliseptol Foam Pure
- ▶ Terralin liquid
- ▶ Ultrazol active

CAUTION**Dirt reduces the light strength**

- ▶ Keep cover clear through regular cleaning.
- ▶ Only wipe cleaning allowed.

- ▶ Clean the PA clear cover using a non-abrasive cleaning cloth and a suitable cleaning agent

CAUTION

- ▶ Wipe cleaning is only allowed in a horizontal position.

CAUTION

To minimize the risk of disease transmission, applicable health and safety regulations and the requirements of the national hygiene and disinfection authorities must be observed in addition to these instructions.

8.1 Sterilizing the hand grip

- Sterilization must be done according to **ISO 17665-1** (Sterilization of health care products in moist heat).

CAUTION

Damage to the hand grip

- Do not sterilize with hot air
- Package the hand grip in a sterile bag before sterilization.
- The hand grip is designed exclusively for damp sterilization with 3 times fractioned pre-vacuum and saturated steam with the following parameters:

Temperature	134°C
Overpressure	2.0 bar
Dwell time	6 min
Vacuum drying	20 min

- After sterilization, check the hand grip for mechanical integrity.
- Do not use damaged hand grips.

9. SAFETY INSPECTIONS

DANGER

Death hazard from electric shock.

- Unplug from the grid and turn switch to off position.
- The power supply cable must be checked at least once a year for damage.

CAUTION

- Maintenance and repairs can be performed only by qualified electricians.
- The corresponding user profile is in Chapter 2 Safety instructions.


EVERY YEAR:

- Check the power supply cable for damage and replace if necessary.
- Check for deformations and cracks in the plastic parts.
- Check the load bearing system for distortion or damage.
- Check for loose parts.
- Check the connection of the M8 safety screw and nut, as well as the grub screws on the ceiling bracket

9.1 Setting the spring force


CAUTION

- The ideal spring force is set at the factory.

- As a first measure slightly tighten the joint cover (by a ¼ turn); this to achieve greater stability of the luminaire head, if necessary

Otherwise:

- Remove the service cover of the spring balanced arm (loosen the 2 crosstip screws)
- Place a suitable tool (e.g. a 3-type Allen key) into the hole 1, turn by ¼ turn in the direction of the arrow, a repeat a number of times until the desired stability is reached
- Screw the service cover on again (attention: do not tighten too strongly)

10. DISMANTLING

DANGER

Death hazard from electric shock.

Before dismantling the luminaire, disconnect all pins from the power supply.

WARNING

Injury hazard

The spring balanced arm is under a high spring force. If the terminal device is not in the uppermost spring balanced arm position, the spring balanced arm will move rapidly upward and cause severe injury. Only dismantle the terminal device when the spring balanced arm is in the uppermost position

10.1 Disposal

Do not discard the luminaire with household waste. Follow local regulations and bring the luminaire to a disposal site or give it to a dealer with appropriate service.

Cut the cable directly at the casing.


The products listed above are over 95% recyclable. For a high percentage of the used materials to be either physically reused or used for energy after the end of this product's life cycle, the luminaires have been designed with recycling in mind. They do not contain materials that are hazardous or require monitoring.


► Disposable cover (order no D15.445.000)

11. ACCESSORIES

► Handle for TRIANGO 130, 100 and 60 (order no D10.442.000)

► Sterilizable surgical handle for TRIANGO 130, 100 and 60 (order no D10.052.000)

12. ADDITIONAL INFORMATION

The luminaire itself is maintenance free.

Additional documents may be requested from the manufacturer for this product.

Using this luminaire does not pose a risk to other equipment. To save energy, the luminaire should be switched on only when it is actually needed.

Any serious incident that has occurred with the product **must be reported** to the manufacturer or their representative and the responsible authorities of the country in which the user is located.

13. TROUBLESHOOTING

Fault	Possible cause	Troubleshooting	User profile
The luminaire does not light up	Contact fault	Switch on again	All
The luminaire does not light up	Lamp defective	Contact manufacturer's service dept.	Only by manufacturer's service dept.
The luminaire does not light up	No mains voltage	Check voltage, check all connections	Qualified Electrician

14. TECHNICAL DATA

Electrical Data:	
Rated input voltage	100 – 240 VAC
Frequency range	50 / 60 Hz
Power consumption	max. 90 VA
Input current	760 – 320 mA
Power factor	0.69 – 0.45
Integral electronic transformer	24 VDC output
Photometric values:	
Central illuminance E_c at 1.0m distance	100'000 lx **
Light field diameter d10 at a distance of 1.0 m: TRIANGO 100-1, TRIANGO 100-3, TRIANGO Endo 100-1 TRIANGO Fokus 100-1, TRIANGO Fokus 100-3	dia. = 18 cm * dia. = 18 cm / 23 cm / 28 cm*
Light field diameter d50 at a distance of 1.0 m:	dia. = 10 cm *
Color temperature: TRIANGO 100-1, TRIANGO Fokus 100-1, TRIANGO Endo 100-1 TRIANGO 100-3, TRIANGO Fokus 100-3	4700K * 3700K*/ 4300K*/ 4700K*
Color rendering Index Ra TRIANGO 100-1, TRIANGO Fokus 100-1, TRIANGO Endo 100-1 TRIANGO 100-3, TRIANGO Fokus 100-3	94* 96*
Color rendering Index R9: TRIANGO 100-1, TRIANGO Fokus 100-1, TRIANGO Endo 100-1 TRIANGO 100-3, TRIANGO Fokus 100-3	72* 92*
Depth of illumination L1 + L2	74 cm *
Relationship between irradiance E_e and irradiance E_c :	3.6 mW/m ² /lx
Remaining illuminance when dimmed by a mask TRIANGO 100-1 Versions TRIANGO Fokus 100-3 Versions	31.6 % 3700 K: 30.7% 4300 K: 26.8% 4700 K: 23.2%
Remaining illuminance when dimmed by two masks TRIANGO 100-1 Versions TRIANGO Fokus 100-3 Versions	49.4 % 3700 K: 56.0% 4300 K: 54.3% 4700 K: 52.5%
Remaining illuminance at the bottom of a simulated cavity TRIANGO 100-1 Versions, TRIANGO Fokus 100-3 Versions	100 %
Remaining illuminance at the bottom of a simulated cavity with one mask TRIANGO 100-1 Versions TRIANGO Fokus 100-3 Versions	31.3% 3700 K: 30.7% 4300 K: 26.8% 4700 K: 23.3%
Remaining illuminance at the bottom of a simulated cavity with two masks TRIANGO 100-1 Versions TRIANGO Fokus 100-3 Versions	49.6% 3700 K: 54.1% 4300 K: 52.0% 4700 K: 50.2%

	TRIANGO 100-1 Versions	TRIANGO Fokus 100-3 Versions
Maximum illuminance distance (DMI)	860 mm	815 mm
Maximum irradiance and radiance settings	Distance 860 mm Maximum illuminance	Distance 815/815/815 mm Maximum illuminance Narrow focus position
Total irradiance	487,5 W/m ²	3700 K: 402,2 W/m ² 4300 K: 446,3 W/m ² 4700 K: 430,9 W/m ²
ratio d50/d10	0.59	0.51 0.55 0.51
		* -10% / +20% tolerance ** only for 4700K / D10 = 180mm
Ambient conditions for transport, storage and operation:		
Ambient temperature (storage and transport)		-20°C to +70°C
Ambient temperature (operation)		+10°C to +35°C
Rel. humidity (non-condensing) (transport, storage and operation)		max. 75%
Maximum use height (operating)		
- standard version		3000m (above sea level)
- (H) version		5000m (above sea level)
Atmospheric pressure (operating)		
- standard version		70-106 kPa
- (H) version		54-106 kPa
Maximum height (storage and transport)		Unlimited
Weight:		
Luminaire head		2.4 kg
TRIANGO 100 C		15.9 kg
TRIANGO 100 W		15.3 kg
TRIANGO 100 F		20.0 kg
Operating mode:		
Operating mode		Continuous operation
Classification:		
TRIANGO 100		Protection class I
Degree of protection according to IEC 60529		IP 20
Luminaire head		IP 43 (horizontal position)
Classification according to EU-REGULATION 2017/745 (MDR), Article 51		Class I
U.S. FDA Device Class		Class I
Electrical safety testing and EMC according to:		AAMI ES60601-1: 2005/A2:2021 CAN/CSA-C22.2 No 60601-1:14 IEC 60601-1:2005 + A2:2020 IEC 60601-1-2:2014 + A2:2020 IEC 60601-2-41:2021
Blue light danger according to EN/IEC 62471		RG 2 (moderate risk)
Life cycle of the light source:		
Life cycle		50,000 h (L70/B50)

15. ELECTROMAGNETIC COMPATIBILITY (EMC)

Medical electrical devices are subject to special precautionary measures with regard to electromagnetic compatibility. This device can be affected by other electrical devices.

WARNING

Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the TRIANGO 100, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

Electromagnetic environment

The device may only be operated in the environments specified in the “Intended Use” section of the Operating Instructions. The medical device is intended for use in the electromagnetic environment specified below

The above listed models are intended for use in the electromagnetic environment specified below. The customer or the user of the above listed models should assure that they are used in such an environment.		
Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The above listed models use RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The above listed models are suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations / flicker emissions IEC 61000-3-3	Complies	

The above listed models are intended for use in the electromagnetic environment specified below.
The customer or the user of the above listed models should assure that they are used in such an environment.

Enclosure Port

Immunity test	Test Condition	IEC 60601 Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV Contact $\pm 2,4,8,15$ kV Air	± 8 kV Contact $\pm 2,4,8,15$ kV Air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Radiated RF EM fields and Proximity fields from RF wireless communications equipment IEC 61000-4-3	3 V/m 80 MHz – 2,7 GHz 80% AM 1kHz	3 V/m	Mains power quality should be that of a typical commercial or hospital environment.
	385MHz (18Hz Pulse Modulation)	27 V/m	
	450MHz (FM+/-5KHz deviation 1kHz sine or 18Hz Pulse Modulation)	28 V/m	
	710MHz (217Hz PM)	9 V/m	
	745MHz (217Hz PM)	9 V/m	
	780MHz (217Hz PM)	9 V/m	
	810MHz (18Hz PM)	28 V/m	
	870MHz (18Hz PM)	28 V/m	
	930MHz (18Hz PM)	28 V/m	
	1720MHz (217Hz PM)	28 V/m	
	1845MHz (217Hz PM)	28 V/m	
	1970MHz (217Hz PM)	28 V/m	
	2450MHz (217Hz PM)	28 V/m	
	5240MHz (217Hz PM)	9 V/m	
	5500MHz (217Hz PM)	9 V/m	
	5785MHz (217Hz PM)	9 V/m	
RATED power frequency magnetic fields IEC 61000-4-8	30 A/m 50Hz or 60Hz	30 A/m	Power frequency magnetic fields should be at levels characteristic of a magnetic field typical location in a typical commercial or hospital environment.
IMMUNITY to proximity magnetic fields	30 kHz CW Test level: 8 A/m 134,2 kHz PM 2,1 kHz Test level: 65 A/m 13,56 MHz PM 50 kHz Test level: 7,5 A/m	Not applicable.	The TRIANGO 100 light do not contain magnetically sensitive components or circuitry within the ENCLOSURE.

The above listed models are intended for use in the electromagnetic environment specified below.
The customer or the user of the above listed models should assure that they are used in such an environment.

Input a.c. power PORT

Immunity test	Test Condition	IEC 60601 Compliance level	Electromagnetic environment - guidance
Electrical fast transient/bursts IEC 61000-4-4	± 2 kV 100kHz Repetition frequency	± 2 kV	Mains power quality should be that of a typical commercial or hospital environment.
Surges IEC 61000-4-5	$\pm 0,5$ kV, ± 1 kV line(s) to line(s) $\pm 0,5$ kV, ± 1 kV, ± 2 kV line(s) to ground	± 1 kV, Differential mode ± 2 kV Common mode	Mains power quality should be that of a typical commercial or hospital environment.
Conducted RF induced by RF fields IEC 61000-4-6	3 V 0,15 MHz – 80 MHz 6 V m) in ISM and amateur radio bands between 0,15 MHz and 80 MHz (n) 80 % AM at 1 kHz	3 V 0,15 MHz – 80 MHz 6 V m) in ISM radio bands between 0,15 MHz and 80 MHz (n) 80 % AM at 1 kHz	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0% U_T ; 0°, 45°, 90°, 135°, 180°, 225°, 270°, 315° 0% U_T ; 0° 0% U_T ; 70% 0% U_T ; 0%	0,5 Cycles 1Cycle 25/30 Cycles (50/60Hz) 250/300 Cycles (50/60Hz) (5s)	Mains power quality should be that of a typical commercial or hospital environment. If the user of the above listed models requires continued operation during power mains interruptions, it is recommended that the above listed models are powered from an uninterruptible power supply or battery.

Comment: n) The ISM (industrial, scientific and medical) bands between 0,15 MHz and 80 MHz are 6,765 MHz to 6,795 MHz; 13,553 MHz to 13,567 MHz; 26,957 MHz to 27,283 MHz; and 40,66 MHz to 40,70 MHz. The amateur radio bands between 0,15 MHz and 80 MHz are 1,8 MHz to 2,0 MHz, 3,5 MHz to 4,0 MHz, 5,3 MHz to 5,4 MHz, 7 MHz to 7,3 MHz, 10,1 MHz to 10,15 MHz, 14 MHz to 14,2 MHz, 18,07 MHz to 18,17 MHz, 21,0 MHz to 21,4 MHz, 24,89 MHz to 24,99 MHz, 28,0 MHz to 29,7 MHz and 50,0 MHz to 54,0 MHz.

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USA

Waldmann Lighting Company
9, W. Century Drive
Wheeling, Illinois 60090
USA
Telephone +1 847 520 1060
Fax +1 847 520 1730
www.waldmannlighting.com
waldmann@waldmannlighting.com

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
Further distribution partners you find at: www.derungs.swiss.

Recommended Installation

TRIANGO FOCUS 130-3

TRIANGO 100-1 C


TRIANGO 100-3 C

TRIANGO FOCUS 100-1


TRIANGO FOCUS 100-3

TRIANGO ENDO 100-1


Installation Height

Installation Height (with intermediate ceiling adapter)

Installation Position

* If installation height = 195cm,
and bed height = 80cm
and lighting distance = 100cm (100'000 lux)

**If possible, we recommend that the light is
installed opposite the treating doctor.
This makes operating the light more ergonomic

Ceiling base plate – Hole pattern