


Portable endoscope

GT-28

- ✓ Portable, ergonomically designed handle
- ✓ Can be connected to mobile phones, tablets, computers
- ✓ Multiple specifications available to meet various diagnosis and treatment needs




Convenient for carrying or onsite service



- 1. lens
- 2. LED
- 3. channel outlet

Portable endoscope



Model	Insertion Tube Diameter	Field of View	Distal End Diameter	Observation Range	Instrument Channel	Bending Capability	Image Output
GT-58(B)	Φ5.8mm	140°	2.6mm	3-2000	650mm	180° up and down	USB to computer or mobile phone etc