



# Compressor Nebulizer AirForce One

The AirForce One compressor nebulizer is a cost effective and reliable device for aerosol therapy at home or in hospitals. Compact, lightweight and quiet in operation, the AirForce One allows discreet aerosol therapy. Engineered with the patient in mind, the high performance yet easy-to-use features offers the user an effective and simple drug delivery option.

## Main Features

- Only weighs 1.5kg AirForce One is the lightest aerosol product in the Drive DeVilbiss International range
- Durable design for everyday use including a 3 year warranty
- Built in carrying handle for added convenience
- Offering a maximum flow rate of 5.1l/min
- Ready to Use! Comes complete with accessories and carry bag

## TECHNICAL SPECIFICATIONS

Dimensions	19.5cm L x 9.5cm H x 15cm W
Weight	1.5kg
Voltage	230VAC, 50HZ
Power	100W
Consumption	
Noise Level	52 dBA
Max Flow	5.1 l/min
Max Pressure	2,2 bar

The AirForce One is available with a UK or Euro power plug.

## UK version comes with a Disposable nebulizer kit

Containing an instruction manual, disposable Evo nebulizer, nose piece, mouthpiece with exhalation valve, air tubing, adult and paediatric mask and 5 x filters.



Technical specifications of Evo nebulizer	
Minimum filling volume	2ml
Maximum filling volume	8ml
Nebulization Rate	0,4ml/min
Particle Size (MMAD)	3,16 µm

## Euro version comes with a Reusable nebulizer kit

Containing an instruction manual, reusable breath enhanced nebulizer kit "MedelJet Plus", nose piece, mouthpiece with exhalation valve, air tubing, adult and paediatric mask and 5 x filters.



Technical specifications of MedelJet Plus	
Minimum filling volume	2ml
Maximum filling volume	8ml
Nebulization Rate	0,3ml/min
Particle Size (MMAD)	< 5 µm

## Drive DeVilbiss International

DeVilbiss Healthcare GmbH  
Kamenzer Straße 3  
68309 Mannheim, Germany

UK +44 (0) 121 / 521 31 40  
DE +49 (0) 621 / 1 78 98 0  
kontakt@devilbisshc.com

[www.drivedevilbiss-int.com](http://www.drivedevilbiss-int.com)

