



Item number 1999915542



Underpad, ABENA Tissue Underpads, 6 layers, 60x40cm, light blue

- ✓ Absorbent core with multilayer tissue
- ✓ Soft and smooth nonwoven surface
- ✓ Rapid and effective liquid absorption
- ✓ Liquid safe back sheet and sealed edges for extra protection
- ✓ Available in different sizes and differing numbers of tissue layers



Product Description

ABENA Tissue Underpads are used as protective disposable pads for e.g. beds and chairs to provide user peace of mind in any situation. The disposable pads are primarily intended for non-sterile procedures such as wound care and blood sampling, but are also suitable for use as protective pads for incontinence problems. The underpads have sealed edges for optimum leakage protection, and the surface is covered with a soft nonwoven material for optimum user comfort. In addition to helping to provide peace of mind when sitting or lying down, the disposable pads also help reduce the risk of contamination and increase hygiene. ABENA Tissue Underpads are available in different sizes and with either six or eight layers of tissue.



Specifications

Base name	Underpad
Brand	ABENA
Sub-brand	Tissue Underpads
Sub-color	Light blue
Features	Non-sterile
Single or multiple use	Single use
Material	Nonwoven, PE, tissue
Length/depth	60 cm
Width	40 cm
Weight, net	45 g
Tolerance	+/- 1 cm
Layer quantity	6 layers
Absorption Rothwell	300 ml
Certifications	CE, MD
CE Class (Medical Devices)	Class I
Directives, regulations and acts	MDR (EU) 2017/745
Because of the products composition it has an expiry date from the production date on:	5 years
Storage Instructions	Store dry, room temperature and no direct sunlight.
Product Disposal Instructions	Can be disposed of with normal household waste sorted according to local regulations. If contamination is suspected, the product must be discarded as clinical waste.
Packaging Disposal Instructions	To be recycled or incinerated.

Instructions for use/application

The underpad must be placed on the bed or chair with the blue side down. Ensure that the underpad is fully smoothed out before use.

Packaging data

Unit	Contains	Length	Width	Height	EAN
cII	16 pck	42.8 cm	36.6 cm	40.8 cm	5713571006948
pck	25 pcs	24.7 cm	24 cm	10.2 cm	5713571006931



Regulation (EU) 2017/745 is a regulation of the European Union on the clinical investigation and sale of medical devices for human use. It repeals Directive 93/42/EEC, which concerns medical devices, and Directive 90/385/EEC, which concerns active implantable medical devices, on 26 May 2021.



The CE mark guarantees that a product is safe to use and complies with all safety precautions. CE stands for Conformité Européenne (European Conformity) and is mainly found on electronic equipment, safety equipment, construction products, and medical devices.



The product is only for single-use.

