

PRODUCT SPECIFICATION

Art. No.	
NAME OF PRODUCT	Pasteur pipette
DESCRIPTION OF PRODUCT	Pasteur pipette , transparent, 1ml, 150mm length
INTENDED USE OF PRODUCT	Used in laboratory

RAW MATERIAL		
CODE FOR THE ENCLOSED	RAW MATERIAL	ADDITIONAL INFORMATION
RAW MATERIAL SPECIFICATION	Pipette	PE
	box	Single Wall Corrugated Fibreboard
	Transport carton	Double Wall Corrugated Fibreboard
	Packing strap	PP+PET

QUALITY SPECIFICATION OF THE PRODUCT			
CHARACTERISTIC	TARGET VALUE	QUALITY TOLERANCES	TEST METHOD
FINISHED SIZE	Total Length=150mm	+/- 0.5mm	Measuring

DESCRIPTION OF THE PACKAGING SYSTEM						
PACKING	BOX SIZE (cm)			CARTON SIZE (cm)		
	L	W	H	L	W	H
500pc//box, 300pcs/ctn	26	16	21	50	27	43

