



RUDAMED[®]-sensitive

Product Description and Purpose

The white first-aid dressing is made of a crosswise elastic non-woven fabric and is provided with a hypoallergenic polyacrylate adhesive. The dressing is especially suited for sensitive skin. RUDAMED[®]-sensitive is packed in a convenient dispenser box with opening.

Composition

Polyethylene, polyester and polyacrylate adhesive

Normative and Legal Requirements

RUDAMED[®]-sensitive is a medical product in accordance with the MDD, class I, rule 4, and complies with the requirements of the directive 93/42/EEC.

The product does not contain dangerous toxic substances according to REACH. It has CE marking and DIN EN ISO 15223-1 labels on all its packaging.

Packaging

Primary packaging: folding box made of cellulose

Secondary packaging: carton made of cellulose

Storage

To be stored in a dry and dust-free environment, protected against direct sunlight

