



## NOBAFILM®

REF 173027

### Product Description and Purpose

The sterile incision film is made of a transparent, hypoallergenic polyurethane film which is coated with a hypoallergenic adhesive. NOBAFILM® is waterproof, but steam and air permeable. The dimensions are 30 cm x 27 cm total size, 24 cm x 27 cm adhesive surface. Microorganisms cannot permeate the film. It is used to cover the operation area in a sterile way. Owing to the fingerlifts on both sides, the surgical site can be optimally covered. The mechanical barrier function minimizes the introduction of germs into the surgical site.

### Contraindications

Known allergy against a component listed below

### Composition

Polyurethane film, polyacrylate adhesive

### Normative und gesetzliche Anforderungen

NOBAFILM® is a medical product according to the MPG, the directive 93/42/EEC and the regulation MDR (EU 2017/745) and is classified as class Is, according to rule 4.

Sterilization of the product complies with DIN EN 11135.

The product does not contain dangerous toxic substances according to REACH.

It has CE marking and DIN EN ISO 15223-1 labels on all its packaging.

### Packaging

Primary packaging:	paper-film packaging
Secondary packaging:	folding box made of cellulose
Tertiary packaging:	carton made of cellulose

### Storage

To be stored in a dry and dust-free environment

The product bears following labels and symbols:



STERILE EO



CE  
0123